

Easy as 1-2-3 instructions on how to install or replace the Cleanviro waterless urinal tablet

1. FOR NEW INSTALLATION, TURN OFF THE MAIN WATER SUPPLY OF THE URINAL AND DRAIN THE WATER TANK BY FLUSHING THE URINAL UNTIL IT IS EMPTY.
2. Using screw driver, remove the stainless strainer from the wall urinal if it is still in place (for new installation).



3. Prepare a bucket of water and put some liquid detergent



4. Pour the water in the urinal drain. The water should settle down the drain quickly. If it doesn't then the urinal drain could be clogged, in that case use a small drain rubber pump (see picture below) to unblock the clogging. Repeat the process until the water settles down quickly – in some cases it takes about 10 buckets to clear up the drain block.



5. Wipe the urinal face using a soft sponge or a damp cloth



6. Pour three capfuls of the cleanviro liquid concentrate (undiluted) onto the urinal drain. Remember to pour concentrated liquid and not a diluted one as it will not be as effective in eliminating old odour buildup in the drain.



7. Put the plastic strainer at the bottom of the urinal



8. Place the cleanviro urinal tablet on top of the plastic strainer



9. Place the ceramic cover to hide the urinal tablet from view



10. Finally spray about 10 times the urinal face and the top of the ceramic cover with the Cleanviro liquid spray cleaning solution. For instructions on how to prepare the liquid cleaning solution, please refer to the label in the Liquid Concentrate bottle.



11. Turn off the water supply. You now have successfully converted your wall urinal into a waterless urinal system. Congratulations.



WHEN TO REPLACE THE URINAL TABLET.

Normally the Cleanviro urinal tablet will last between 2-4 weeks depending on use and conditions (heavier usage means shorter efficacy). You can inspect the tablet once a week and if you notice it's already significantly smaller than the drain hole, or when it's losing its effectiveness (urine smell becoming noticeable) then it's time to replace the tablet.

Easy as 1-2-3 instructions on daily maintenance of urinal with Cleanviro waterless urinal tablet

A. THROUGH OUT THE DAY

1. IT IS STRONGLY ADVISED TO SPRAY THE URINAL WITH THE CLEANVIRO LIQUID SPRAY CLEANING SOLUTION TWICE DAILY TO MAINTAIN THE EFFECTIVENESS OF THE WATERLESS URINAL.
2. Spray the urinal face and the top of the ceramic cover with the Cleanviro liquid spray cleaning solution. Ideally spray the urinal and its surroundings about 15 times. You can also spray the wall and floor to maintain hygienic atmosphere.



3. Wipe the urinal face using a soft sponge or a damp cloth



4. Ideally you do this simple task once during the heavy traffic period (during or just after lunch time) and once in the afternoon, but you can spray the solution as often as you want depending on the condition.



B. END OF THE DAY MAINTENANCE

1. Using gloves, remove the ceramic cover, plastic strainer and existing tablet and quickly rinse it with water.



2. Prepare a bucket of water and put some liquid detergent



3. Pour the water in the urinal drain.



4. Wipe the urinal face using a soft sponge or a damp cloth



5. Pour three capfuls of the cleanviro liquid concentrate (undiluted) onto the urinal drain. Remember to pour concentrated liquid and not a diluted one as it will not be as effective in eliminating old odour buildup in the drain.



6. Put the plastic strainer at the bottom of the urinal



7. Place the used cleanviro urinal tablet on top of the plastic strainer



8. Place the ceramic cover to hide the urinal tablet from view



9. Finally spray about 10 times the urinal face and the top of the ceramic cover with the Cleanviro liquid spray cleaning solution. For instructions on how to prepare the liquid cleaning solution, please refer to the label in the Liquid Concentrate bottle.



That's it – you just completed the daily maintenance.